

SKIN STAPLER

Patent

No.ZL201310353428.X



• **Safety Lock**
Patent No.ZL200920233463.7



REMOVER

SKIN STAPLER

Patent

No.ZL201310353428.X



Safety Lock

Patent No.ZL200920233463.7

SKIN STAPLER

Features:

- 1、Staples are compatible with skin tissue well, no rejection reaction;
- 2、Support different dimensions' staples, apply to stapling different skin wounds;
- 3、Much faster than suturing by hand, thus relieve operation workload, shorten operation time and decrease occurring of complication;
- 4、Small stapled wounds, avoid injuring skin tissue again, thus benefit wounds healing;
- 5、Staple beautifully, smoothly and accurately, healing scar is small and neat, when remove staples, no painful feeling;
- 6、Angled head of stapler provides clear view to ensure precise stapling.

Types:

Unit: mm

Model	Type	Pin (staple)			Quantity Of Pins	Optional (Pin remover)
		Diameter	Length	Height		
KHPF	W1	0.58	13.05	3.5	10、 15、 25、 35、 45	KHPF- CW
	W2					KHPF- CW
	R	11.4	2.7	KHPF- CB		

JIANGSU COOPWIN MED S&T CO., LTD

Address : B3 Hutang Science & Technology Industrial Park, NO198, Xiahe Road, Hutang Town, Wujin Area, Changzhou City, P.R.C

Tel: +86-0519-68212728 Fax : +86-0519-68212727

Website : www.coopwinmed.com

Email: sale1@coopwinmed.com ; sale2@coopwinmed.com ; feedback@coopwinmed.com

TROCAR



- **No-conversion universal sealing cap**

Support 5mm-12mm diameters' operation instruments

- **Specific thread make sure excellent stability**

- **Blade/bladeless/ optical types support**

TROCAR

Features:

- 1、 Excellent air tightness of trocar, well protect pneumoperitoneum during operation;
- 2、 The cannula surface has specific thread, make sure excellent stability, thus reduce trauma to abdominal wall;
- 3、 The cannula enable 3mm-15mm diameters' operation instruments to go through;
- 4、 No need replacing switch cap during operation;
- 5、 More easier to control during operation,one hand using only.

KHTC-Blade types:

Unit:mm

Models	Inner diameter of pipe, D	Length of pipe, L
	Dimension	Dimension
KHTC-5/65	5.9	64
KHTC-5/100	5.9	97
KHTC-10/100	10.5	98
KHTC-12/100	12.9	98
KHTC-15/100	15.6	110

KHCC-Bladeless types:

Unit:mm

Models	Inner diameter of pipe, D	Length of pipe, L
	Dimension	Dimension
KHCC-5/100	5.9	100
KHCC-10/100	10.9	100
KHCC-12/100	13	100

JIANGSU COOPWIN MED S&T CO., LTD


Address : B3 Hutang Science & Technology Industrial Park, NO198, Xiahe Road, Hutang Town, Wujin Area, Changzhou City, P.R.C

Tel: +86-0519-68212728 Fax : +86-0519-68212727

Website : www.coopwinmed.com

Email:sale1@coopwinmed.com ; sale2@coopwinmed.com ; feedback@coopwinmed.com

BIOPSY FORCEPS — KHTS



Oval cup



Oval cup




Alligator cup



Alligator cup

BIOPSY FORCEPS — KHFB A(B)



Oval cup



Oval cup



Alligator cup



Alligator cup

DISPOSABLE ENDOSCOPIC LINEAR CUTTER STAPLER



DISPOSABLE ENDOSCOPIC LINEAR CUTTER STAPLER

Features:

- 1、 Knife inside the cartridge, Better incision effect and avoid cancer cell transfer
- 2、 Space locking device, Make sure the staplesto be formed consistently and good hemostasis result
- 3、 Safety of empty cartridge, Prevent mis-operation of empty cartridge before firing
- 4、 Much larger hatch, Convenient tissue placement and adjustment
- 5、 Low profile anvil design, Convenient equipment placement, reduce tissue damage
- 6、 6 rows of compact ranged staples, Provide good hemostasia result
- 7、 All reloads could be used in the same instrument,Reduce complicated operations, reduce patient costs

Indication:

Transection, resection, and the creation of anastomoses for digestive tissue and organ in endoscopic surgery.

Types:

Models	Size	Unit:mm	
		L	L2
KHQJ -60		280(±30)	60(±15)
KHQJ -160		380(±30)	160(±15)
KHQJ -260		480(±30)	260(±15)

Reloads:

Models	L1	Unit:mm			
		Endoscopic linear cutter cartridges DER			Tolerance
		L3			
KHQJR-30G	30(±3.0)	2.0	2.0	2.0	±0.2
KHQJR-30T		2.0	2.5	3.0	
KHQJR-30P		3.0	3.5	4.0	
KHQJR-45G	45 (±3.0)	2.0	2.0	2.0	±0.2
KHQJR-45T		2.0	2.5	3.0	
KHQJR-45P		3.0	3.5	4.0	
KHQJR-45B		4.0	4.5	5.0	
KHQJR-60T	60(±3.0)	2.0	2.5	3.0	±0.2
KHQJR-60P		3.0	3.5	4.0	
KHQJR-60B		4.0	4.5	5.0	

JIANGSU COOPWIN MED S&T CO., LTD


Address : B3 Hutang Science & Technology Industrial Park, NO198, Xiahe Road, Hutang Town, Wujin Area, Changzhou City, P.R.C

Tel: +86-0519-68212728 Fax : +86-0519-68212727

Website : www.coopwinmed.com

Email:sale1@coopwinmed.com ; sale2@coopwinmed.com ; feedback@coopwinmed.com

DISPOSABLE LINEAR CUTTER STAPLER AND RELOAD



DISPOSABLE LINEAR CUTTER STAPLER AND RELOAD

Features:

- 1、 Knife inside the cartridge, Better incision effect and avoid cancer cell transfer
- 2、 Space locking device, Make sure the staples to be formed consistently and good hemostasis result
- 3、 Safety of empty cartridge, Prevent mis-operation of empty cartridge before firing
- 4、 Two ways Push- knob, Convenient for operation by left or right hand
- 5、 Special rear lock , realizes one- handed operation, convenient for operation
- 6、 Quick release knob, Unlock the instrument quickly
- 7、 Larger front fork hatch, convenient for operation
- 8、 The proximal preloading device eliminates the difference in scissors and makes the nail height uniform
- 9、 Ergonomic anti-slip design for improved safety, stability and comfort in handling the instrument

Indication:

Transection, resection, and the creation of anastomoses for digestive tissue and organ

Types:

Models	Specification	L1		Unit:mm
		Basic size	Tolerance	
KHLC	60 (3.8、 4.6)	62	±2	
	80 (3.8、 4.6)	82		
	100 (3.8、 4.6)	101		

Reloads:

Models	Specification	Nail number	Height	Tolerance	Unit:mm
KHLRC	60 (3.8、 4.6)	64	3.8 4.6	±0.3	
	80 (3.8、 4.6)	84			
	100 (3.8、 4.6)	104			

JIANGSU COOPWIN MED S&T CO., LTD


Address : B3 Hutang Science & Technology Industrial Park, NO198, Xiahe Road, Hutang Town, Wujin Area, Changzhou City, P.R.C

Tel: +86-0519-68212728 Fax : +86-0519-68212727

Website : www.coopwinmed.com

Email:sale1@coopwinmed.com ; sale2@coopwinmed.com ; feedback@coopwinmed.com

DISPOSABLE LINEAR STAPLER AND RELOAD



DISPOSABLE LINEAR STAPLER AND RELOAD

Features:

- 1、 Much larger parallel jaw, Contains more tissue
- 2、 Two segment close, Convenient to adjust the tissue
- 3、 Auto/Manual to push the tissue locating pin, Convenient for operation
- 4、 Full specifications, More options for digestive system operations
- 5、 Low profile anvil design, Convenient equipment placement, reduce tissue damage
- 6、 Human body engineering design, Convenient to operate by one hand
- 7、 Automatically clamps tissue, making one-handed firing easier
- 8、 Imported titanium nails for higher strength and stretch
- 9、 The height of the staple cartridge can be adjusted according to the thickness of the tissue

Indication:

Creation of anastomoses, suture the end and incision for digestive tissue and organ

Types:

Unit:mm

Models	Size	L1	Tolerance	H	Tolerance	Models	Size	L1	Tolerance	H	Tolerance			
KHZX A-30×3.8		30	±2.0	3.8	±2.0	KHZX B-30×3.8		30	±2.0	3.8	±2.0			
KHZX A-45×3.8		45				4.8	45	KHZX B-45×3.8					45	4.8
KHZX A-60×3.8		60						KHZX B-60×3.8					60	
KHZX A-90×3.8		91						KHZX B-90×3.8					91	
KHZX A-30×4.8		30		KHZX B-30×4.8						30				
KHZX A-45×4.8		45		KHZX B-45×4.8			45							
KHZX A-60×4.8		60		KHZX B-60×4.8			60							
KHZX A-90×4.8		91		KHZX B-90×4.8			91							

Reloads:

Unit:mm

Models	Size	L1	Tolerance	H	Tolerance	Models	Size	L1	Tolerance	H	Tolerance
KHZXR A-30×3.8		30	±2.0	3.8	±2.0	KHZXR B-30×3.8		30	±2.0	3.8	±2.0
KHZXR A-30×4.8		30		4.8		KHZXR B-30×4.8		30		4.8	
KHZXR A-45×3.8		45		3.8		KHZXR B-45×3.8		45		3.8	
KHZXR A-45×4.8		45		4.8		KHZXR B-45×4.8		45		4.8	
KHZXR A-60×3.8		60		3.8		KHZXR B-60×3.8		60		3.8	
KHZXR A-60×4.8		60		4.8		KHZXR B-60×4.8		60		4.8	
KHZXR A-90×3.8		90		3.8		KHZXR B-90×3.8		90		3.8	
KHZXR A-90×4.8		90		4.8		KHZXR B-90×4.8		90		4.8	

JIANGSU COOPWIN MED S&T CO., LTD

Address : B3 Hutang Science & Technology Industrial Park, NO198, Xiahe Road, Hutang Town, Wujin Area, Changzhou City, P.R.C

Tel: +86-0519-68212728 Fax : +86-0519-68212727

Website : www.coopwinmed.com

Email:sale1@coopwinmed.com ; sale2@coopwinmed.com ; feedback@coopwinmed.com

DISPOSABLE CIRCULAR STAPLER



DISPOSABLE CIRCULAR STAPLER

Features:

- 1、Integrated Anvil, Make sure better “ B” shaped staple
- 2、Larger inner diameter, Reduce the risk of anastomosis straitness
- 3、Special design for stomach lead hole, Convenient for operation in esophagus
- 4、Imported titanium staples, Make sure stable anastomosis
- 5、More compact staple range, Provide good hemostasia

Indication:

End to end, end to side, side to side anastomosis for esophagus, stomach and intestines

Types:

Size		Unit:mm							
Models	L	Tolerance	D	Tolerance	D1	Tolerance	H	Tolerance	
KHGXA-17×4.2	455	±30	17	±1.0	10	±1.0	4.2	±0.3	
KHGXA-21×4.5			21		13		4.5		
KHGXA-23×4.5			23		15		4.5		
KHGXA-26×4.5			26		17		4.5		
KHGXA-29×4.8			29		21		4.8		
KHGXA-32×4.8			32		22		4.8		
KHGXA-34×4.8			34		24.5		4.8		
KHGXB-17×4.0			17		10		4.0		
KHGXB-21×4.8			21		13		4.8		
KHGXB-21×4.0			21		13		4.0		
KHGXB-23×4.8			23		14.5		4.8		
KHGXB-23×4.0			23		14.5		4.0		
KHGXB-26×4.8			26		16.5		4.8		
KHGXB-26×4.0			26		16.5		4.0		
KHGXB-29×4.8			29		19.5		4.8		
KHGXB-29×4.0			29		19.5		4.0		
KHGXB-32×4.8			32		22.5		4.8		
KHGXB-34×4.8			34		24.5		4.8		

JIANGSU COOPWIN MED S&T CO., LTD

Address : B3 Hutang Science & Technology Industrial Park, NO198, Xiahe Road, Hutang Town, Wujin Area, Changzhou City, P.R.C

Tel: +86-0519-68212728 Fax : +86-0519-68212727

Website : www.coopwinmed.com

Email:sale1@coopwinmed.com ; sale2@coopwinmed.com ; feedback@coopwinmed.com

HEMORRHOIDS STAPLER



Accessory of tpe KHGCF J-M



TST surgical accessory of tpe KHGCF J-L



HEMORRHOIDS STAPLER

Features:

- 1、Integrated Anvil, Make sure better " B" shaped staples
- 2、Integrated shaft, Provide more intensity & stability
- 3、Much bigger & deeper capacity of cartridge, Contains more tissue
- 4、Imported Titanium staples, Provide much higher strength & stretch
- 5、More staples, better hemostasis
- 6、Using technologically innovative materials, the sound is clearer, making the operation safer
- 7、Accurate minimum gap to avoid tissue damage caused by excessive compression

Indication:

For the selective removal of mucosa on the dentate line

Product sterilization method:

Sterilized with ethylene oxide, the product should be sterile and disposable

Types:

Models	O.D	I.D	Unit:mm
			Height
KHGCA-32x3.8	32	22.8	3.8
KHGCA-34x3.8	34	24.8	
KHGCA-36x3.8	36	26.8	
KHGCA-32x3.8	32	22.8	
KHGCA-34x3.8	34	24.8	
KHGCA-36x3.8	36	26.8	

JIANGSU COOPWIN MED S&T CO., LTD


Address : B3 Hutang Science & Technology Industrial Park, NO198, Xiahe Road, Hutang Town, Wujin Area, Changzhou City, P.R.C

Tel: +86-0519-68212728 Fax : +86-0519-68212727

Website : www.coopwinmed.com

Email: sale1@coopwinmed.com ; sale2@coopwinmed.com ; feedback@coopwinmed.com

HEMORRHOIDAL MULTI-BAND LIGATOR WITH ANOSCOPE



Accessory

HEMORRHOIDAL MULTI-BAND LIGATOR WITH ANOSCOPE

Features:

- 1、 Automatic, time-saving, less effort, simple and practical
- 2、 Only 5- 10 minutes are used for one person to complete the operation.
- 3、 Clinic treats without anesthesia.
- 4、 Slightly pains, rare complications.

Application:

- 1、 Various hemorrhoids
- 2、 Mixed hemorrhoids
- 3、 For persons who have hemorrhoids block or incomplete anal pad rebound after PPH or other therapy, this product can be used as a supplementary treatment.
- 4、 Other: rectal focal disease, such as rectal polyps, rectum hemangioma (vascular malformation), etc.

Indication:

Suitable for ligation of various hemorrhoids, mixed hemorrhoids and rectal benign polyps in each period

Types:

Unit:mm

Models	Adsorption cartridge outer diameter	Rubber gasket inner diameter
KHTZA-2/A(B)	12.5	2.2
KHTZA-3/A(B)		
KHTZA-4/A(B)		
KHTZA-5/A(B)		
KHTZB 2/A(B)	12	2
KHTZB-3/A(B)		
KHTZB-4/A(B)		

JIANGSU COOPWIN MED S&T CO., LTD

Address : B3 Hutang Science & Technology Industrial Park, NO198, Xiahe Road, Hutang Town, Wujin Area, Changzhou City, P.R.C

Tel: +86-0519-68212728 Fax : +86-0519-68212727

Website : www.coopwinmed.com

Email:sale1@coopwinmed.com ; sale2@coopwinmed.com ; feedback@coopwinmed.com

DISPOSABLE CIRCUMCISION STAPLER



Standardization , safety and simpleness of the surgery operation

Shorter surgery time (3-5mins)

Automatic suturing , less pain

Mild edema, little bleeding, few complications

Rapid recovery, no need return for second time

Attractive appearance, high satisfaction